

# SHARPS Transport Tube™



**Sharps Transport Tube™**, manufactured by Sharps Compliance, Inc., has been created for simplified management of used sharps, such as needles, syringes, and lancets. Sharps® Transport Tubes are the ideal solution for “on-the-go” individuals who self inject drug therapies and travel, live or work in environments where short-term containment of used sharps is necessary. Transport tubes are not designed to be mailed once full; nor should they be put into the trash. Place them in a sharps container intended for proper disposal, such as the Sharps Disposal by Mail System® container.

The Sharps Disposal by Mail System® is the perfect solution for safe disposal of used needles, lancets, and syringes including those that have been contained in the Sharps Transport Tube. Keeping these potentially infectious sharps out of the trash protects our children, pets, and the environment.

## SHARPS TRANSPORT TUBE

### FEATURES

- ♦ Thick, puncture resistant, two piece, rigid polypropylene material
- ♦ Easy grip shape for transport
- ♦ Ready to use
- ♦ Environmentally safe, may be autoclaved or incinerated
- ♦ Screw-on top ensures safe transport

### SPECIFICATIONS

- Capacity: 1-2 syringes, up to 12 pen needles, or up to 8 lancets
- Product #50035

## HOW THE SHARPS DISPOSAL BY MAIL SYSTEM® WORKS

The **Sharps Disposal By Mail System**® is a government approved system comprised of a specially designed mailable sharps container, 4-mil bag liner, postage pre-paid shipping box, multi-part manifest tracking form, and instructions for use. The cost of mailing and destruction is included in the purchase price.

When the sharps container is full, simply place it back into the pre-paid shipping box. Remove and complete the tracking form located in the plastic pouch on the outside of the shipping box. Keep the bottom copy and return the remaining copies to the pouch. Give the package to your postal carrier. It will be mailed directly to our treatment facility and destroyed without ever opening the box.



Don't forget to maintain your supply of Sharps Transport Tubes and an extra Sharps Disposal by Mail System® so you are never caught without a safe disposal method.

**View Industry-related tutorials now available on our website: [www.sharpsinc.com](http://www.sharpsinc.com).**

### How the System Works

1. Order the Sharps Disposal By Mail System® online or by phone. If you have a business account, be sure to sign up for SharpsTracer®.



2. Deposit sharps and other medical waste into appropriate container. Close lid securely when ready for disposal.



3. Complete manifest tracking form and keep Generator copy.



4. Package container in the Sharps postage pre-paid return mailback box following enclosed instructions.



5. Give sealed box to your mail carrier, take to a U.S. Post Office, or request a Priority Pickup from the Post Office.



6. Medical waste is delivered to the Sharps dedicated treatment facility. Your waste is weighed, documented and completely destroyed without ever opening the box.



7. Proof of destruction is available online by logging in to your SharpsTracer® account.

