




Pitch-It® IV POLE SYSTEM

The Problem

Traditional IV poles are heavy, cumbersome, and can limit patient use and mobility. Not only does the weight and size of these poles increase initial shipping costs, having to retrieve, clean, and prepare them for reuse can be expensive and time consuming as well.



The Pitch-It® IV Pole Solution

Sharps' patented Pitch-It® IV poles are designed to eliminate the frustration of retrieving, cleaning, bagging, and tagging old IV poles. The lightweight Pitch-It® IV poles can be packaged with other patient supplies with significant reduction in weight of the shipment.



The poles are quickly unfolded, and their innovative tripod design allows them to be used on any floor or tabletop. Patients can use them at home or fold them up to take anywhere they need to continue their therapy.



Pitch-It® IV poles also eliminate the need for home infusion patients to improvise gravity-fed infusions with such items as coat hangers, door knobs, light fixtures or even nails in the wall. When no longer needed, the aluminum pole can be set out with other recyclables.



ADVANTAGES

- Lower shipping costs
- Loved by both nurses and patients
- Portable and collapsible
- Lightweight
- Convenient for transport
- Durable
- Easy to assemble
- No cleaning
- Disposable – place in recycling
- No pickups required

SPECIFICATIONS

PITCH-IT® JR.

- Item #30004-012
- Tabletop model for gravity-fed iv infusions
- Extended height 53" (1309 MM)
- Collapsible height: 28"
- Weight: 1 pound

PITCH-IT®

- Item #30003-012
- Floor Model for gravity-fed IV infusions
- Extended height 72" (1779 mm)
- Weight: 1 pound

PITCH-IT® SR.

- Item #30002-006
- Durable enough for most infusion pumps
- Two telescoping sections that twist and lock to adjust the pole at different heights
- Five legs for easy movement and less chance of tipping
- Designed to hold up to 3 liters of IV fluid
- Easily identifiable pump mounting collar
- Extended height: 74" (1829 mm)
- Weight: 4 pounds

SHARPS
Defining Change Through Innovation...
Sharps Disposal By Mail™

9220 Kirby Drive Suite 500
Houston, Texas 77054
800.772.5657 713.432.0300
Fax: 713.432.0555
www.sharpsinc.com