



# SHARPS DISPOSAL BY MAIL SYSTEM®



**Sharps Disposal By Mail System®** is the ideal disposal solution for:

- Syringes discovered at your location
- Sharps and small quantities of other medical waste generated by your facilities
- Syringes used by self-injectors at home or elsewhere



The Sharps Disposal by Mail System® is government-approved and includes a specially designed mailable sharps container, a postage prepaid return shipping box with barcode tracking, a protective 3-mil bag liner, instructions in both English and Spanish, and a four-part simplified manifest tracking form.



The cost of mailing and disposal, as well as confirmation of destruction, is included in the up-front purchase price.

When the Sharps Disposal by Mail System® sharps container is full, it is placed in the pre-paid postage return box and delivered via the U.S. Postal Service to the Sharps treatment facility in Carthage, Texas. Upon receipt, Sharps documents the delivery, weight and destruction of each container and its contents. Confirmation of destruction is made available electronically using Sharps proprietary tracking system SharpsTracer®.



## ADVANTAGES








- Sharps collection, shipment and disposal all in one system and all for one price
- Containers resist punctures and leaks
- No waiting for pickups and risking non-compliance
- Systems are ordered when needed—no contracts or minimum service levels
- No hidden fees or additional charges
- No big box to store
- No interruptions during patient care
- Packaged containers are simply given to a mail carrier for shipment

## WHAT MAKES THE SHARPS DISPOSAL BY MAIL SYSTEM® EVEN BETTER

- In use since 1994
- Sharps does it all from manufacturing through destruction
- One company and one disposal site for all locations
- USPS-permitted and approved by all 50 states
- System complies with OSHA, EPA and all federal regulations
- Nationally recognized team of regulatory experts and compliance consultants at your disposal
- Online manifest tracking and management with SharpsTracer®

Syringe Capacity Per Container			
Item #	Size	3cc syringes	1cc syringes
10100	1-Quart	20-30	50-70
11000	1-Gallon	95-100	210-230
12000	2-Gallon	250-270	440-460
13000	3-Gallon	450-480	700-720
15000	5-Gallon	For disposal of red bag biohazardous waste. <b>Not for disposal of loose sharps.</b> Sharps must be contained in OSHA-compliant sharps containers.	

## HOW THE SHARPS DISPOSAL BY MAIL SYSTEM® WORKS

How the System Works	
1. Order the Sharps Disposal By Mail System® online or by phone. If you have a business account, be sure to sign up for SharpsTracer®.	
2. Deposit sharps and other medical waste into appropriate container. Close lid securely when ready for disposal.	
3. Complete manifest tracking form and keep Generator copy.	
4. Package container in the Sharps postage pre-paid return mailback box following enclosed instructions.	
5. Give sealed box to your mail carrier, take to a U.S. Post Office, or request a Priority Pickup from the Post Office.	
6. Medical waste is delivered to the Sharps dedicated treatment facility. Your waste is weighed, documented and completely destroyed without ever opening the box.	
7. Proof of destruction is available online by logging in to your SharpsTracer® account.	

**How it Works!**  
**View Industry-related video**  
**tutorials available on our website**  
**[www.sharpsinc.com](http://www.sharpsinc.com).**