

1.5-Quart Sharps Recovery System Multipacks

Six 1.5-Quart #11506, Twelve 1.5-Quart #11512

Perfect for offices with multiple treatment locations, Multipack Sharps Recovery Systems include multiple sharps containers with the convenience of a single mailback box.

The Sharps Recovery System is an innovative system designed for simplified and environmentally-friendly treatment of used needles, syringes and other regulated medical waste.

SYSTEM INCLUDES

- Sharps collection containers
- Single prepaid USPS return shipping box
- Waste manifest tracking form
- Proper treatment and online documentation of treated waste

ADVANTAGES

- Cost-effective alternative to pickup service, reducing costs by as much as 50% or more
- No contracts or hidden fees - order systems as needed
- Approved for sharps and other small quantities of medical waste, keeping you in compliance
- Easy to package for return, saving you time for other important job responsibilities
- Option to schedule free USPS priority pickup online at your convenience
- Online tracking and proof of treatment via SharpsTracer available 24/7, giving you complete cradle-to-grave protection



Twelve 1.5-Quart #11512

Six 1.5-Quart Sharps Recovery System Multipack

Product Specifications and Shipping Data

#11506 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	4.5" L	4.5" W	8.625" H
Single Mailback Box:	15.125" L	10.25" W	10" H
Shipping Box:	15.375" L	10.75" W	11" H
UPC Code:	634188000097		
Max Return Weight:	25 lbs		
Units / Pallet:	48		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	1152 cases		

Twelve 1.5-Quart Sharps Recovery System Multipack

Product Specifications and Shipping Data

#11512 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	4.5" L	4.5" W	8.625" H
Single Mailback Box:	19.5" L	14.875" W	10.125" H
Shipping Box:	19.875" L	15.375" W	11" H
UPC Code:	634188000103		
Max Return Weight:	25 lbs		
Units / Pallet:	24		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	576 cases		

Sharps Secure System

#50060

- One 1.5-1uart sharps collection container*
- One stainless steel cabinet
- Two keys

**Please Note: Item #50060 does not include a mailback shipping box. The Sharps Recovery System (#10150), which includes a replacement 1.5-qt. sharps container and a mailback shipping box, must be purchased separately.*



For more information, please call **800.772.5657**
or visit us online at www.sharpsinc.com.

© 2015 Sharps Compliance, Inc. All rights reserved.

★ MADE IN ★
USA