

## 20- & 28-Gallon Sharps Recovery Systems

20-Gallon #10020, 28-Gallon #12800

**For a cost-effective and regulatory-compliant disposal solution for sharps containers and sealed red bag waste, the 20- and 28-Gallon Sharps Recovery Systems are ideal.**

The Sharps Recovery System is an innovative system designed for simplified and environmentally-friendly treatment of regulated medical waste. This leak-proof, puncture-resistant system accommodates your smaller existing sharps containers. It also features a locking pin for when it is ready to return for treatment in order to secure your medical waste during transport.

### SYSTEM INCLUDES

- 20 or 28-gallon collection container
- Prepaid USPS return shipping box
- Waste manifest tracking form
- Proper treatment and online documentation of treated waste

### ADVANTAGES

- Cost-effective alternative to pickup service, reducing costs by as much as 50% or more
- No contracts or hidden fees – order systems as needed
- Approved for sharps containers (no loose sharps) and other small quantities of medical waste, keeping you in compliance
- Easy to package for return, saving you time for other important job responsibilities
- Option to schedule free USPS priority pickup online at your convenience
- Online tracking and proof of treatment via SharpsTracer available 24/7, giving you complete cradle-to-grave protection



20-Gallon #10020

#### 20-GALLON SYSTEM ACCOMMODATES THE FOLLOWING SHARPS CONTAINERS

1-Quart:	1.5-Quart:	5-Quart:	1-Gallon:	2-Gallon:	3-Gallon:
up to 28	up to 18	up to 6	up to 8	up to 5	up to 3

#### 28-GALLON SYSTEM ACCOMMODATES THE FOLLOWING SHARPS CONTAINERS

1-Quart:	1.5-Quart:	5-Quart:	1-Gallon:	2-Gallon:	3-Gallon:
up to 35	up to 28	up to 8	up to 12	up to 7	up to 4

## 20-Gallon Sharps Recovery System

### Product Specifications and Shipping Data

#### #10020 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	21.5" L	21.25" W	18.5" H
Single Mailback Box:	22" L	22" W	20" H
Shipping Box:	22" L	22" W	20" H
UPC Code:	634188100209		
Max Return Weight:	35 lbs		
Units / Pallet:	12 units		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	288 cases		

## 28-Gallon Sharps Recovery System

### Product Specifications and Shipping Data

#### #12800 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	19" L	20" W	23.125" H
Single Mailback Box:	20" L	20" W	24" H
Shipping Box:	20" L	20" W	24" H
UPC Code:	634188128005		
Max Return Weight:	35 lbs		
Units / Pallet:	8 units		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	192 cases		

## Rolling Base

#80005

Easily maneuver larger systems throughout your facility with the steel Rolling Base, perfect for the 20- and 28-Gallon Systems. The rolling base features four heavy-duty casters and locking wheels.



For more information, please call **800.772.5657** or visit us online at [www.sharpsinc.com](http://www.sharpsinc.com).

© 2015 Sharps Compliance, Inc. All rights reserved.

★ MADE IN ★  
**USA**