

20- & 30-Gallon TakeAway Recovery Systems

20-Gallon #80020, 30-Gallon#80030

For a cost-effective and regulatory-compliant disposal solution for sharps containers and used healthcare materials, the 20- and 30-Gallon TakeAway Recovery Systems are ideal.

The TakeAway Recovery System is an innovative system designed for simplified and environmentally-friendly treatment of sharps containers and other used healthcare materials. This leak-proof, puncture-resistant system accommodates your smaller existing sharps containers. It also features a locking pin for when it is ready to return for treatment in order to secure your sharps during transport.

SYSTEM INCLUDES

- Sharps collection container
- Prepaid UPS return shipping box
- Tracking form
- Proper treatment and online documentation of treated materials

ADVANTAGES

- Cost-effective alternative to pickup service, reducing costs by as much as 50% or more
- No contracts or hidden fees - order systems as needed
- Approved for sharps containers (no loose sharps) and other used healthcare materials, keeping you in compliance
- Easy to package for return, saving you time for other important job responsibilities
- Online tracking and proof of treatment via SharpsTracer available 24/7, giving you complete cradle-to-grave protection



30-Gallon#80030

20-GALLON SYSTEM ACCOMMODATES THE FOLLOWING SHARPS CONTAINERS

1-Quart:	1.5-Quart:	5-Quart:	1-Gallon:	2-Gallon:	3-Gallon:
up to 28	up to 18	up to 6	up to 8	up to 5	up to 3

30-GALLON SYSTEM ACCOMMODATES THE FOLLOWING SHARPS CONTAINERS

1-Quart:	1.5-Quart:	5-Quart:	1-Gallon:	2-Gallon:	3-Gallon:
up to 40	up to 32	up to 10	up to 14	up to 8	up to 5

20-Gallon TakeAway Recovery System

Product Specifications and Shipping Data

#80020 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	21.5" L	21.25" W	18.5" H
Single Mailback Box:	22" L	22" W	20" H
Shipping Box:	22" L	22" W	20" H
UPC Code:	634188800208		
Max Return Weight:	50 lbs		
Units / Pallet:	12		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	288 cases		

30-Gallon TakeAway Recovery System

Product Specifications and Shipping Data


#80030 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	20.5" L	20.3" W	28" H
Single Mailback Box:	22" L	22" W	30.75" H
Shipping Box:	22" L	22" W	30.75" H
UPC Code:	634188800307		
Max Return Weight:	50 lbs		
Units / Pallet:	8		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	192 cases		

Rolling Base

#80005

Easily maneuver larger systems throughout your facility with the steel Rolling Base, perfect for the 20- and 30-Gallon Systems. The rolling base features four heavy-duty casters and locking wheels.



For more information, please call **800.772.5657** or visit us online at www.sharpsinc.com.

© 2015 Sharps Compliance, Inc. All rights reserved.