

28-Gallon Preloaded Sharps Recovery System

Ten 1-Gallon #12810, Six 5-Quart #12861, Six 2-Gallon #12862

For a cost-effective and regulatory-compliant disposal solution for sharps containers and sealed red bag waste, the 28-Gallon Preloaded Sharps Recovery Systems are ideal for offices with multiple treatment locations.

The Sharps Recovery System is an innovative system designed for simplified and environmentally-friendly treatment of sharps containers and other regulated medical waste. This large leak-proof, puncture-resistant container comes preloaded with multiple smaller sharps containers.

SYSTEM INCLUDES

- Sharps collection containers
- Single prepaid USPS return shipping box
- Waste manifest tracking form
- Proper treatment and online documentation of treated waste

ADVANTAGES

- Cost-effective alternative to pickup service, reducing costs by as much as 50% or more
- No contracts or hidden fees - order systems as needed
- Approved for sharps containers and tied red bag waste, keeping you in compliance
- Easy to package for return, saving you time for other important job responsibilities
- Option to schedule free USPS priority pickup online at your convenience
- Online tracking and proof of treatment via SharpsTracer available 24/7, giving you complete cradle-to-grave protection



Six 2-Gallon #12862

28-Gallon Medical Professional Sharps Recovery System with Ten 1-Gallon Sharps Collection Containers

Product Specifications and Shipping Data

#12810 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	19" L	20" W	23.125" H
Single Mailback Box:	20" L	20" W	24" H
Shipping Box:	20" L	20" W	24" H
UPC Code:	634188000714		
Max Return Weight:	35 lbs		
Units / Pallet:	8		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	192 cases		



Wall Mount Brackets for 1-, 2- and 3-Gallon Sharps Collection Containers

#50004, #50005, #50014

Conveniently mount 1-, 2- and 3-gallon sharps containers with three bracket styles: basket, front and side.

28-Gallon Medical Professional Sharps Recovery System with Six 5-Quart Sharps Collection Containers (Red)

Product Specifications and Shipping Data

#12861 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	19" L	20" W	23.125" H
Single Mailback Box:	20" L	20" W	24" H
Shipping Box:	20" L	20" W	24" H
UPC Code:	634188000134		
Max Return Weight:	35 lbs		
Units / Pallet:	8		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	192 cases		

28-Gallon Medical Professional Sharps Recovery System with Six 2-Gallon Sharps Collection Containers

Product Specifications and Shipping Data

#12862 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	19" L	20" W	23.125" H
Single Mailback Box:	20" L	20" W	24" H
Shipping Box:	20" L	20" W	24" H
UPC Code:	634188000707		
Max Return Weight:	35 lbs		
Units / Pallet:	8		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	192 cases		



Rolling Base

#80005

Easily maneuver larger systems throughout your facility with the steel Rolling Base, perfect for the 20- and 28-Gallon Systems. The rolling base features four heavy-duty casters and locking wheels.

For more information, please call **800.772.5657** or visit us online at www.sharpsinc.com.

© 2015 Sharps Compliance, Inc. All rights reserved.