

## 30-Gallon Preloaded TakeAway Recovery System

Ten 1-Gallon #83101, Six 5-Quart #83061, Six 2-Gallon #83062

**For a cost-effective and regulatory-compliant disposal solution for sharps containers and sealed red bag waste, the 30-Gallon Preloaded TakeAway Recovery Systems are ideal for offices with multiple treatment locations.**

The TakeAway Recovery System is an innovative system designed for simplified and environmentally-friendly treatment of sharps containers and other used healthcare materials. This large leak-proof, puncture-resistant container comes preloaded with multiple smaller sharps containers. It also features a locking pin for when it is ready to return for treatment in order to secure your used healthcare materials during transport.

### SYSTEM INCLUDES

- Sharps collection containers
- Single prepaid UPS return shipping box
- Tracking form
- Proper treatment and online documentation of treated materials

### ADVANTAGES

- Cost-effective alternative to pickup service, reducing costs by as much as 50% or more
- No contracts or hidden fees - order systems as needed
- Approved for sharps containers and other small quantities of used healthcare materials, keeping you in compliance
- Easy to package for return, saving you time for other important job responsibilities
- Option to schedule free UPS priority pickup online at your convenience
- Online tracking and proof of treatment via SharpsTracer available 24/7, giving you complete cradle-to-grave protection



Ten 1-Gallon #83101

## 30-Gallon TakeAway Recovery System with Ten 1-Gallon Sharps Containers

Product Specifications and Shipping Data

#83101 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	20.5" L	20.33" W	28" H
Single Mailback Box:	22" L	22" W	30.75" H
Shipping Box:	22" L	22" W	30.75" H
UPC Code:	634188831011		
Max Return Weight:	50 lbs		
Units / Pallet:	8		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	192 cases		



## Wall Mount Brackets for 1-, 2- and 3-Gallon Sharps Collection Containers

#50004, #50005, #50014

Conveniently mount 1-, 2- and 3-gallon sharps containers with three bracket styles: basket, front and side.

## 30-Gallon TakeAway Recovery System with Six Covidien 5-Quart Sharps Containers

Product Specifications and Shipping Data

#83061 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	20.5" L	20.33" W	28" H
Single Mailback Box:	22" L	22" W	30.75" H
Shipping Box:	22" L	22" W	30.75" H
UPC Code:	634188830618		
Max Return Weight:	50 lbs		
Units / Pallet:	8		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	192 cases		

## 30-Gallon TakeAway Recovery System with Six 2-Gallon Sharps Containers

Product Specifications and Shipping Data

#83062 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	20.5" L	20.33" W	28" H
Single Mailback Box:	22" L	22" W	30.75" H
Shipping Box:	22" L	22" W	30.75" H
UPC Code:	634188830625		
Max Return Weight:	50 lbs		
Units / Pallet:	8		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	192 cases		



## Rolling Base

#80005

Easily maneuver larger systems throughout your facility with the steel Rolling Base, perfect for the 20- and 30-Gallon Systems. The rolling base features four heavy-duty casters and locking wheels.

For more information, please call **800.772.5657**  
or visit us online at [www.sharpsinc.com](http://www.sharpsinc.com).

© 2015 Sharps Compliance, Inc. All rights reserved.

★ MADE IN ★  
**USA**