

## 5-Quart Sharp Recovery System

5-Quart #10501, 5-Quart Intro Kit #10501-INTRO

**Small quantity generators in need of a cost-effective and regulatory-compliant disposal solution for sharps will find the 5-Quart Sharps Recovery System the perfect fit for their needs.**

The 5-quart sharps containers' automatic shut off prevents overfilling, while the bold red "FULL" indicates when to change the container. The anti-kickback lid limits accidental or intentional access to its contents, and the locking tabs provide extra safety during final closure. The wall enclosure is ideal for securing sharps containers, keeping counter space clear and helping prevent access to contents in high traffic areas.

### SHARPS RECOVERY SYSTEM (#10501) INCLUDES

- One Covidien® 5-quart sharps container
- One prepaid USPS return shipping box
- Waste manifest tracking form
- Proper treatment and online documentation of treated waste

### INTRO KIT (#10501-INTRO) INCLUDES

- Two Covidien® 5-quart sharps containers
- One wall enclosure/wall mount unit
- One prepaid USPS return shipping box
- Waste manifest tracking form
- Proper treatment and online documentation of treated waste

### ADVANTAGES

- Cost-effective alternative to pickup service, reducing costs by as much as 50% or more
- No contracts or hidden fees - order systems as needed
- Approved for sharps and other small quantities of medical waste, keeping you in compliance
- Easy to package for return, saving you time for other important job responsibilities
- Option to schedule free USPS priority pickup online at your convenience
- Online tracking and proof of treatment via SharpsTracer available 24/7, giving you complete cradle-to-grave protection

*COVIDIEN, COVIDIEN with logo, and the Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. © Covidien. All Rights Reserved*



5-Quart #10501

**5-Quart Sharps Recovery System**

Product Specifications and Shipping Data

#10501 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	12.25" L	4.75" W	10.5" H
Single Mailback Box:	14" L	7.75" W	15.25" H
Shipping Box:	14.75" L	8.25" W	16" H
UPC Code:	634188105013		
Max Return Weight:	10 lbs		
Units / Pallet:	48		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	1152 cases		

**5-Quart Sharps Recovery System Intro Kit**

Product Specifications and Shipping Data

#10501-INTRO PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	12.25" L	4.75" W	10.5" H
Single Mailback Box:	14" L	7.75" W	15.25" H
Shipping Box:	14.25" L	13.5" W	16" H
UPC Code:	634188105006		
Max Return Weight:	10 lbs		
Units / Pallet:	27		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	648 cases		

**Wall Enclosure System\***

#10502

- One Covidien® 5-quart sharps container
- One wall enclosure/wall mount unit

\*Please Note: System does not include a shipback box.



For more information, please call **800.772.5657**  
or visit us online at [www.sharpsinc.com](http://www.sharpsinc.com).