

Sharps Accessories

Scoop Cones #50002-012, Finger Cones #50001-014, Transport Tubes #50035-006

Scoop Cones

Scoop Cones offer a hands-free pickup option for sharps. They have a flat surface on one edge of the opening with an attached lid for scooping needles directly into the cone. Cones are designed for one-time use and can be properly disposed of in Sharps Recovery Systems.

Finger Cones

Finger Cones are similar to Scoop Cones but have a round opening with no flat edge. Needles are placed into cones and then disposed of in larger sharps containers.

Transport Tubes

Transport Tubes have been created for simplified management of used sharps, such as needles, syringes, and lancets. They are the ideal solution for “on-the-go” individuals who self-inject drug therapies and travel, live or work in environments where short-term containment of used sharps is necessary. Transport tubes are not designed to be mailed or put into the trash. Place them in a sharps container intended for proper disposal, such as the Sharps Recovery System.



Finger Cones #50001-014

Scoop Cones

Product Specifications and Shipping Data

#50002-012 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	8.5" L	2.5" W	2.25" H
Shipping Box:	19.25" L	4.25" W	2.5" H
UPC Code:	63418800059		
Units / Pallet:	418		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	10032 cases		

Finger Cones

Product Specifications and Shipping Data

#50001-014 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	8.5" L	2.5" W	2.25" H
Shipping Box:	18" L	11.75" W	18.25" H
Units / Pallet:	16		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	384 cases		

Transport Tubes

Product Specifications and Shipping Data

#50035-006 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	2.25" L	8.25" W	1.75" H
Shipping Box:	9" L	11.75" W	2.75" H
UPC Code:	634188500351		
Units / Pallet:	256		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	6144 cases		