

Two 3.5-Gallon Pail TakeAway Recovery System

Two 3.5-Gallon Pails #83502

Replace expensive waste hauling services with the TakeAway Recovery Pail, a cost-effective, UPS shipback system.

The TakeAway Recovery System is an innovative system designed for simplified and environmentally-friendly treatment of used needles, syringes and other used healthcare materials.

SYSTEM INCLUDES

- Two 3.5-gallon pails
- One prepaid UPS return shipping box
- Waste manifest tracking form
- Proper treatment and online documentation of treated waste

ADVANTAGES

- Cost-effective alternative to pickup service, reducing costs by as much as 50% or more
- No contracts or hidden fees – order systems as needed
- Approved for sharps and other used healthcare materials, keeping you in compliance
- Easy to package for return, saving you time for other important job responsibilities
- Online tracking and proof of treatment via SharpsTracer available 24/7, giving you complete cradle-to-grave protection



Two 3.5-Gallon Pails #83502



Two 3.5-Gallon Pail TakeAway Recovery System

Product Specifications and Shipping Data

#83502 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	14" L	14" W	10.75" H
Single Mailback Box:	16" L	16" W	22.25" H
Shipping Box:	16.5" L	16.25" W	23" H
UPC Code:	634188835026		
Max Return Weight:	50 lbs		
Units / Pallet:	20		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	480 cases		

For more information, please call **800.772.5657**
or visit us online at **www.sharpsinc.com**.

© 2015 Sharps Compliance, Inc. All rights reserved.

