

## 3-Gallon Sharps Recovery System Multipacks

Two 3-Gallon #13002, Four 3-Gallon #13004

**Perfect for offices with multiple treatment locations, Multipack Sharps Recovery Systems include multiple sharps containers with the convenience of a single mailback box.**

The Sharps Recovery System is an innovative system designed for simplified and environmentally-friendly treatment of used needles, syringes and other regulated medical waste.

### SYSTEM INCLUDES

- Sharps collection containers
- Single prepaid USPS return shipping box
- Waste manifest tracking form
- Proper treatment and online documentation of treated waste

### ADVANTAGES

- Cost-effective alternative to pickup service, reducing costs by as much as 50% or more
- No contracts or hidden fees - order systems as needed
- Approved for sharps and other small quantities of medical waste, keeping you in compliance
- Easy to package for return, saving you time for other important job responsibilities
- Option to schedule free USPS priority pickup online at your convenience
- Online tracking and proof of treatment via SharpsTracer available 24/7, giving you complete cradle-to-grave protection



Two 3-Gallon #13002

**Two 3-Gallon Sharps Recovery System Multipack**

Product Specifications and Shipping Data

#13002 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	8.75" L	5.5" W	16.75" H
Single Mailback Box:	9.5" L	12.5" W	18" H
Shipping Box:	10.25" L	13" W	20" H
UPC Code:	634188130022		
Max Return Weight:	25 lbs		
Units / Pallet:	48		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	1152 cases		

**Four 3-Gallon Sharps Recovery System Multipack**

Product Specifications and Shipping Data

#13004 PRODUCT SPECIFICATIONS / SHIPPING DATA

Single Container:	8.75" L	5.5" W	16.75" H
Single Mailback Box:	18.625" L	12.5" W	18" H
Shipping Box:	19.375" L	13.75" W	19.625" H
UPC Code:	634188130046		
Max Return Weight:	25 lbs		
Units / Pallet:	21		
Total / Pallet / TL:	24 pallets		
Number of Cases / TL:	504 cases		

Wall Mount Brackets for 1-, 2-  
and 3-Gallon Sharps  
Collection Containers

#50004, #50005, #50014

Conveniently mount 1-, 2- and 3-gallon sharps containers with three bracket styles: basket, front and side.



For more information, please call **800.772.5657**  
or visit us online at **www.sharpsinc.com**.

© 2015 Sharps Compliance, Inc. All rights reserved.

★ MADE IN ★  
**USA**